

EntericBio

Certificate

Product Number **EB003**
Lot Number **EB90509**
Expiry Date **2010-05-28**
Date of Issue **2009-08-15**

Optical & Physical Properties

| Test | Specification | Result | Pass / Fail |
|------------|---------------|---------|-------------|
| Appearance | Clear | Clear | Pass |
| | 2 Phase | 2 Phase | Pass |

Quality Test*

| O.D. culture, dilution | Specification (PCR result) | Result | Pass / Fail |
|------------------------|----------------------------|--------|-------------|
| O.D. 0,3-0,4, 1:10 | + | + | Pass |
| O.D. 0,3-0,4, 1:100 | + | + | Pass |
| O.D. 0,3-0,4, 1:1000 | - or + | +/- | Pass |

**Quality test performed acc. to SOP QS-QA-2131-xx (Quality control of EntericBio Extract) using E. coli ATCC 35218 and a PCR protocol of Hsu and Tsen (2001).*

Reference:

Hsu, S.-C. and Tsen, H.-Y. (2001): PCR primers designed from malic acid dehydrogenase gene and their use for detection of Escherichia coli in water and milk samples. Int. J. Food Microbiol. 64: 1-11

Tom Russell.