

EntericBio

Extract

EntericBio Extract removes inhibitors from faecal samples,
To give a simple, fast & effective DNA extraction for PCR applications

Removes PCR Inhibitors

- such as bile salts, cholic acid, deoxycholic acids from faecal samples.

Fast & Effective DNA extraction

- one step – heating and sonication
- heating inactivates nucleases, safeguarding DNA
- sonication ensures fast and effective process

One Step Extraction

- simply add sample and sonicate
- inhibitors concentrated in top phase

Non Toxic Reagents

- chemical free lysis method

More effective than standard methods

- simplified extraction process
- more efficient extraction for *Shigella* & *Campylobacter*

Add Enriched Faecal Sample



Single tube technology, one test four results



Serosep Limited

Unit F5, Eastway Business Park
Ballysimon Road, Limerick
Ireland

Tel (+353) 61 440207 Fax (+353) 61 440237
Email info@serosep.com