

# EntericBio

## Extract

EntericBio Extract removes inhibitors from faecal samples,  
To give a simple, fast & effective DNA extraction for PCR applications

### Removes PCR Inhibitors

- such as bile salts, cholic acid, deoxycholic acids from faecal samples.

### Fast & Effective DNA extraction

- one step – heating and sonication
- heating inactivates nucleases, safeguarding DNA
- sonication ensures fast and effective process

### One Step Extraction

- simply add sample and sonicate
- inhibitors concentrated in top phase

### Non Toxic Reagents

- chemical free lysis method

### More effective than standard methods

- simplified extraction process
- more efficient extraction for *Shigella* & *Campylobacter*

Add Enriched Faecal Sample



*Single tube technology, one test four results*



**Serosep Limited**

Unit F5, Eastway Business Park  
Ballysimon Road, Limerick  
Ireland

Tel (+353) 61 440207 Fax (+353) 61 440237  
Email [info@serosep.com](mailto:info@serosep.com)